





## KR 1000 Infant Radiant Warmer









NOVOS offers a combined warming bed. KR 1000 which is designed for different needs of a newborn in delivery room, emergency and NICU. This minimizes the number of intervention and reduces the stress as a result. Our KR 1000 model infant warmer comes with three operation modes to create the best thermal support, so that newborn can use their energy for growth rather than temperature stabilization. Optional modules and features provide caregiver easy access and utmost capability for a complete support of newborn fortheir different needs.

#### **Advanced Features**

- Easy to use control panel with large and bright LED displays.
- Pre-Warm, Manual and Baby (Servo) operation modes for different purpose of use.
- Controlled heating of infant bed provides a warm welcome to the precious newborn by pre-warm mode.
- SmartSilence helps clinician to keep their hand hygienic and to prevent stress on newborn during any surgical intervention or medical examination by hands free silencing of audible alarms.
- 360° Rotating bed is a comfortable option for both the newborn and clinicians with the possibility of easily rotating the bed platform instead of changing their positions or repositioning the newborn.
- SmartTilt and SmartMemory let clinicians to apply trendelenburg and reverse trendelenburg positions just by pushing of a button.
- Warmer head hosts the most efficient heating system with ceramic heaters.
- Pivoting heater module provides less intervention thus less stress for newborn during X-Ray procedures.
- Digital Appar timer of 1, 5, 10, 20 minutes for appar scoring.
- · Intelligent alarms for maximum patient safety.
- Large mattress and foldable side panels for easy access to newborn.
- Examination lamps for easy observation of newborn.
- Optional measurement of ambient temperature lets caregivers to obtain the ideal warming therapy.

### Optimal Heating and Maximum Patient Safety

KR 1000 infant warmer is equipped with an intelligent control panel with microprocessor for continuous optimal thermal support to newborn and easy usage for caregiver. Perfectly designed heater module guarantees an equal distribution of heat over the mattress securing equal heat distribution throughout newborn body. For optimal thermoregulation of newborn, KR 1000 measures newborn skin temperature with a dual thermistor skin probe. During operation, newborn skin temperature measured by dual thermistor skin probe to maximize safety and accuracy of newborn warming. For maximum patient safety, system checks itself and simultaneously informs the user by auto-pausing thermo regulation and intelligent safety precaution alarms.







#### **USER FRIENDLY DESIGN**

KR 1000 offers you the comfort and freedom that you need, while you give the best care to your precious patients.

#### **Control Panel**

Easy to use control panel provides ability to make easy adjustments of the parameters. You can easily read all the necessary information on the large, bright LED displays.

#### **Examination Lamps**

Two angled examination lamps provides a clear observation of newborn.

#### SmartSilence and Intelligent Alarm System

Smart silence feature provides caregiver hand-free silencing of the alarm thuskeeping clinician's hands hygienic during operation. Intelligent alarm system warns caregiver if an unexpectedly long operation time is passed without reaching target skin temperature in baby mode. Also in manual mode a similar intelligent alarm is triggered if the device is operated over safety limits.

#### Pivoting Heater Module

Pivoting heater module allows easy integration of a mobile x-ray unit over the mattress, which means less intervention and less stress for newborn. Also thanks to pivoting heater module, doctors or clinicians can reach the patient for examination purpose easily.

#### Foldable Side Panels

Easy patient access from three sides of the bed is ensured with foldable side panels.

#### X-Ray Cassette Tray

X-ray casette tray is accessible from front side of the bed.

#### SmartTilt and SmartMemory Systems

Thanks to the SmartTilt feature, trendelenburg and reverse trendelenburg positioning of KR 1000 bed is handled by an electronical control system with a simple push of a button located on the front side of the bed unit. SmartMemory System offers two memory slots available for positioning the bed in two different angles with one touch easily.



















#### **FULLY CONFIGURABLE FOR MULTIPLE NEEDS**


KR 1000 is designed to fulfill all requirements and needs with the following optional modules.



Optional 360° Rotating bed provides easy access to newborn during surgical operations. With this feature clinicians can easily rotate the bed platform instead of changing their positions or repositioning the newborn without any stress.



It is possible to adjust height of the bed module to get a more ergonomical setup with an easy to use foot controlled electrical elevation system.



For newborn who needs respiratory support, Easypuff resuscitation module enables you to resuscitate newborn in a simple and safe way eliminating the risk of barotrauma or insufficient ventilation. You can precisely adjust the PIP and PEEP values, also ensure the ideal resuscitation is obtained.

#### Air & Oxygen Flowmeter System

You can deliver Air & Oxygen to the newborn with the desired oxygen concentration using air & oxygen flowmeter system. This feature ensures protection of newborn from the hazards of oxygen surplus.

#### Suction Module

Unique design of KR 1000 can be combined with two suction alternatives such as regular vacuum driven suction and medical oxygen driven venturi suction. Negative pressure applied can be easily adjusted and continuously monitored.

#### SmartWeigh System

SmartWeigh system guarantees precise weight measurement of the newborn at any point of infant bed. No more need for an external scale.

























#### **Technical Specifications**

**Dimensions** 

Width 90cm Depth 113cm Height (standard) 177cm

Height (adjustable, optional) 177 cm 197 cm (max.)

Mattress Size 70 cm x 55 cm

Maximum Capacity of bed 10 kg Mattress Height (standard) 93 cm

Mattress Height 93 cm-113 cm (max.)

(with opt. bed module elevation)

Characteristics

**Examination Light** 

Supply Voltage 220 VAC (±10%) Supply Frequency 50-60 Hz ±12° **Bed Tilting** LED 2x5W

Set Temperature (Baby Mode) 34–38°C in 0.1°C increments

0-100% in 10% increments Heater Output

Skin Temperature Display 18°C-42°C Temperature Measurement Accuracy ± 0.1°C Hands Free Alarm Silence 10 minutes Pre-Warm Mode w/ Air Probe Option 40% ≤ 25°C

30% > 25°C

Pre-warm Mode Standard Device 40% constant

Easypuff T-Piece Resuscitation Module (Optional)

Peak Inspiratory Pressure (PIP) 5-70 cmH2O @ 8 lpm Positive-end Expiratory Pressure (PEEP) 1-10 cmH2O @ 8 lpm

Adjustable Max. Pressure Limit (preset) 40 cm H2O Flow Capacity 5 lpm-15 lpm Airway Pressure Manometer -20 to 80 cmH2O Display Range (±2% of full scale)

Air & O2 Blender (Optional)

21-100% 02 O2 Concentration Range ±3% full scale O2 Accuracy Flowmeter Range 0-15 lpm

Air & Oxygen Flowmeter (Optional)

Wall Supply Pressure 3-5 bars Cylinder Pressure 3-5 bars Patient Gas Supply Flow Range 0-15 lpm

360° Rotating Bed (Optional)

Suction Module (Optional)

-200 to 0 mmHg Manometer Range

500 ml Safety Jar Capacity

O2 for Positive Pressure Driven Input Gas Type

Venturi Type Suction

Vacuum for Negative Pressure

Driven Regular Suction

Weighing Scale (Optional)

Weighing Range 0-10 kg Weight Display Resolution 1 g Accuracy ±4 g

Integrated LED Phototherapy (Optional)

Intensity  $> 40 \mu \text{W/cm}^2/\text{nm}$ Lamp Usage 50.000 hours

Alarms

High Skin Temperature Low Skin Temperature High Temperature Control Baby

Sensor Disconnection/Failure

Axis Failure System Failure Power Failure

**Order List** 

KR 1000 Infant Radiant Warmer SmartWeigh

Integrated Scale System Height Adjustment

Bililed LED Phototherapy System with Flexible Arm

Easypuff Integrated T-piece Resuscitator

360° Rotating Bed Venturi Type Suction Regular Type Suction Low Flow Blender

Air Oxygen Flowmeter System

Humidifier Bottle Disposable Skin Probe Easy Sliding Drawer Monitor Tray IV Pole

Cylinder Rack (E Type)

Air Probe



Proven quality and 100 % customer satisfaction

Having a very strict quality policy brought us 100 percent customer satisfaction. Novos will continue to design and manufacture high quality products for newborn care; all you have to do is just focusing on your patients.

NOVOS Medical reserves the right to make changes without notice in design, specifications and models

The quality management system at NOVOS Medical Systems is certified according to ISO 13485 and product is certified in accordance with Medical Device Directive (93/42/EEC)

Novos Tıbbi Cihazlar Sanayi Ticaret İthalat Ve İhracat Limited Şirketi İvedikosb Mahallesi 1518. Cad. Matbaacılar Sitesi Sit. No: 2/39 Yenimahalle, Ankara, Turkey

Tel: +90 312 384 15 88 Fax: +90 312 384 15 98